

Product Overview

Description:	Adult Arm Sling
Type:	Adjustable
Colours:	Black
Structure:	Breathable neoprene construction with hook-and-loop (Velcro) closures for adjustable support and comfort.



Description

Provides secure, comfortable, and adjustable support for the left or right arm, elbow, and wrist during recovery from fractures, sprains, or post-surgical procedures. Lightweight yet durable, it offers excellent immobilisation while ensuring all-day comfort for the patient.

Indications

- Stabilises the shoulder, arm and wrist after injury or surgery
- For fractures, rotator cuff injuries, sprains and dislocations

Features

- Breathable design for airflow and reduced heat build-up
- Adjustable shoulder strap for a quick, personalised fit
- Removable webbing swathe strap for added immobilisation
- Soft foam neck/shoulder strap for enhanced comfort
- Ambidextrous design—fits left or right arm
- Hook-and-loop closures for easy on/off
- Lightweight design for day-long wear
- Machine washable to keep clean and fresh

Sizes

	Sling Length	Max Sling Depth	Neck Strap Length	Waist Strap Length
Small	35cm	18cm	90cm	105cm
Medium	40cm	18cm	100cm	110cm
Large	45cm	18cm	125cm	125cm
XLarge	50cm	18cm	130cm	135cm
2XLarge	55cm	18cm	145cm	140cm

Important Notice and Disclaimer

This product is intended to provide support and assistance to individuals with mobility impairments. Before using the product, we strongly recommend that you consult with a qualified healthcare professional to ensure that the device is appropriate for your individual needs and medical condition. It is also important to regularly check the device for any signs of wear or damage.

It is important to use this device with caution and to take appropriate safety measures, such as using appropriate footwear and ensuring that the device is properly adjusted and maintained.

By purchasing and using the assistive mobility device/s listed on this product insert page, you acknowledge that you have read and understood this disclaimer and assume all responsibility for any risks associated with the use of the device.

For more information visit us online or contact us.

