

Product Overview

Description:	Junior and Adult Pneumatic Walker Boot
Type:	Adjustable
Colours:	Black with Grey Sole
Structure:	High-impact plastic shell with soft, breathable foam and fabric lining, secured by durable nylon webbing and hook-and-loop (velcro) straps.



Sizes

	Women's UK	Men's UK
Small	4 - 6	2 - 6
Medium	5.5 - 9.5	5.5 - 9.5
Large	8.5 - 11.5	8.5 - 11.5
XLarge	11.5+	11.5+

Description

Designed to provide superior stabilisation, comfort, and protection for individuals recovering from left or right foot or ankle injuries. Featuring an advanced pneumatic (air) system, this boot offers customisable compression to aid healing, reduce swelling, and enhance patient comfort.

Indications

- Conservative management and post-surgical rehabilitation of Achilles tendon rupture
- Stable fractures of the foot, ankle or lower leg
- Severe ankle sprains or ligament injuries
- Post-operative support for foot, ankle or calf procedures
- Use during trauma recovery and rehabilitation phases

Features

- Anti-slip sole provides enhanced safety and stability
- Lightweight construction enhances patient comfort and ease of movement
- Integrated pneumatic air system inflates via a built-in pump for personalised compression and secure fit
- Air release valve positioned below the pump on either side for easy pressure adjustment
- Soft, breathable inner lining ensures comfort during extended wear
- Removable front plastic guard for added protection; includes detachable foam toe padding and two optional square inserts for customised cushioning
- High-cut design offers full circumferential support and immobilisation
- Shock-absorbing sole minimises heel impact and enhances walking comfort
- Ambidextrous design—fits left or right foot

Important Notice and Disclaimer

This product is intended to provide support and assistance to individuals with mobility impairments. Before using the product, we strongly recommend that you consult with a qualified healthcare professional to ensure that the device is appropriate for your individual needs and medical condition. It is also important to regularly check the device for any signs of wear or damage.

It is important to use this device with caution and to take appropriate safety measures, such as using appropriate footwear and ensuring that the device is properly adjusted and maintained.

By purchasing and using the assistive mobility device/s listed on this product insert page, you acknowledge that you have read and understood this disclaimer and assume all responsibility for any risks associated with the use of the device.

For more information visit us online or contact us.



www.astrogroup.co.za

Leading Aluminium

SUPPLIER | FABRICATOR | SURFACE FINISHER



Ver 2 (16 Apr 2026)