

S-Bend Walking Stick





Product Overview

Description: Walking stick, comfortable for medium to long term use

Handheld, adjustable Type:

Tube finish: Anodised matt 15 microns (μ) Colours: Available in natural colour only

Plastics:

Structure: Adjustable aluminium tube

Leg adjustment: Adjustable to 10 different lengths (25mm apart)

Punch options: Single; or double by arrangement

Classification: ISO 9999

Code: 120306



General Information

| Grip Angle | 90° |
|---------------------|-------------|
| Mass | 485.5g |
| Max. Weight of User | 100kg |
| Leg Section Length | 720 - 940mm |

Fittings

| Handle | Made of polypropylene | |
|--------|--|--|
| | Fitted with moulded super soft grip | |
| Foot | Design for maximum adherence | |
| | Fitted inside with anti-perforation ring | |
| | Replaceable | |

Upper Support (Aluminium 6063-T6)

| Length | | 525mm |
|-------------------|--|------------|
| Outside Diameter | | 25.4mm |
| Mass | | 160g |
| Height Adjustment | | X 10 holes |
| Plastics | Fitted with plastic slide at lower end for noise and vibration reduction | |

Lower Tube (Aluminium 6063-T6)

| Foot Length | | 470mm |
|-------------------|---|--------------------------------|
| Outside Diameter | | 21.6mm |
| Mass | | 163g |
| Height Adjustment | | Fitted with plated spring clip |
| Plastics | Fitted with plug at upper end for noise and vibration reduction | |
| Foot | Fitted with 22mm foot plug at lower end | |

Important Notice and Disclaimer

This product is intended to provide support and assistance to individuals with mobility impairments. Before using the product, we strongly recommend that you consult with a qualified healthcare professional to ensure that the device is appropriate for your individual needs and medical condition. It is also important to regularly check the device for any signs of wear or damage.

It is important to use this device with caution and to take appropriate safety measures, such as using appropriate footwear and ensuring that the device is properly adjusted and maintained.

By purchasing and using the assistive mobility device/s listed on this product insert page, you acknowledge that you have read and understood this disclaimer and assume all responsibility for

For more information visit us online or contact us.



